Only for in vitro use in clinical laboratory

CLINICAL SIGNIFICANCE

Store at 2-8°C

Homocysteine (Hcy) is a thiol-containing amino acid produced by the intracellular demethylation of methionine. Total homocysteine (tHcy) represents the sum of oxidised, protein bound and free forms of Hcy.

Elevated levels of tHcv have emerged as an important risk factor in the assessment of cardiovascular disease¹⁻³. Excess Hcy in the blood stream may cause injury to arterial vessels due to its irritant nature, and result in inflammation and plaque formation, which may eventually cause blockage of blood flow to the heart. Elevated levels of tHcy are also linked with Alzheimers disease⁴ and osteoporosis⁵.

PRINCIPLE OF THE TEST:

Oxidised homocysteine is reduced to free homocysteine. Free Hcy is converted to cvstathionine by the use of CBS (cvstathionine beta-synthase) and excess serine. The cystathionine is then broken down to homocysteine, pyruvate and ammonia. Pyruvate is converted to lactate via lactate dehvdrogenase with NADH as coenzyme. The rate of NADH conversion to NAD+ (ΔA_{340nm}) is directly proportional to the concentration of homocysteine.

WARNINGS AND PRECAUTIONS:

- When handling samples, take great care to prevent infection by HBV, HIV, and HCV, wear rubber gloves.

- Be sure not to allow foreign substances including dust, fungi, bacteria, and detergent to get mixed into samples. In addition, take care to prevent contamination of reagents and cuvettes.

- On testing, wear disposable gloves and avoid oral pipetting to prevent infection

- Dispose of the reagent under a large amount of running water, as each component reagent contains sodium azide. If they get into the eyes or mouth or adhere to skin, first aid measures such as thorough flushing with water should be taken. Consult a physician if necessary.

INSTRUMENTS:

Instrument applications are available upon request.

COMPONENT COMPOSITION:

Component	Ingredients	Concentration in Tests
Reagent 1	Tris Buffer	
_	LDH	35KU/L
	L-Serine	0.76 mmol/L
	TCEP	0.5 mmol/L
	NADH	0.47 mmol/L
	PRESERVATIVES	
Reagent 2	Tris Buffer	
	Cystathionine β-Synthase	20KU/L
	Cystathionine β-lyase	10KU/L
	L-Serinelyase	20KU/L
	PRESERVATIVES	

REAGENT PREPARATION AND STABILITY:

Reagent is ready for use.

If stored at 2-8°C and handled properly, component is stable until expiry date stated on the label

TYPE OF SPECIMEN:

Use serum or heparin plasma as specimen.

It is recommended to follow NCCLS procedures (or similar standardised conditions) regarding specimen handling. Specimen should be collected in an appropriate sample container, with proper specimen identification.

 Serum/plasma should be separated from cells within 8 hours after collection. Stability: 2 weeks at 2-8°C.

TEST PROCEDURE:

Materials required but not supplied HCY Controls and Calibrator General Laboratory Equipment

Assay procedure:

Wavelength:	340 nm
Temperature:	30°C or 37°C
Optical path:	1 cm light path.

STANDARD PROCEDURE:	Blank	Calibrator	Sample			
Reagent 1	960 µl	960 µl	960 µl			
Sample			52 µl			
Calibrator		52 µl				
Gently mix and Incubate at 37°C for 5 minutes						
Reagent 2	260 µl	260 µl	260 µl			
Gently mix and Incubate at 37°C for 1 minutes, then measure the change of						

Optical Density per minute (AOD/min) over the next 2 minutes

Calculation:

Concentration = $\Delta Abs/min_{sample}$ x Concentration of Calibrator

 $\Delta Abs / min_{Calibrator}$

Quality Control

All clinical laboratories should establish an internal Quality Control program. Verify instrument and reagent performance with recommended controls or similar. The values obtained for QC should fall within manufacturer's acceptable ranges or should be established according to the laboratory's QC program.

EXPECTED VALUES:

Adult: ≤15µmol/L Elder population ≥ 60 years: 15 - 20µmol/L

Each laboratory should establish its own reference range. Amylase results should always be reviewed with the patient's medical examination and history.

PERFORMANCE CHARACTERISTICS:

Performance results can vary with the instrument used. Data obtained in each individual laboratory may differ from these values.

Linearity:

This assay is linear up to 50µmol/L. For samples with a higher concentration, dilute 1:1 with 0.9% NaCl (9g/l) and re-assay. Multiply result by 2.

Interfering substances: 5

Results of study are as follows:

Bilirubin (mixed isomers):	Less than 10% interference up to 600 µmol/l Bilirubin							
Haemolysis: Haemoglobin	Less	than	10%	interference	up	to	500	mg/dl
Lipemia:	Less than 10% interference up to 500 mg/dl Lipemia						mia	

Sensitivity:

The Lower Detectable Level was estimated at 0.7µmol/L.

Precision:							
Within Run N = 20	Mean <i>(</i> µmol/L <i>)</i>	SD	% CV	Between Run N = 20	Mean <i>(</i> µmol/L <i>)</i>	SD	% CV
Level 1 Level 2	12.2 25.6	1.02	2.62 1.78	Level 1 Level 2	12.9 26.4	0.96 3.84	2.68 1.92
20.012	23.0	2.72		20.012	23.4	0.04	

Method Comparison:

Using 23 samples, a comparison, between this HCY test (y) and another commercially available test (x), gave the following results:

y = 0.97x - 3.67 r = 0.997 Sample range: 3 to 36µmol/L

BIBLIOGRAPHY:

- 1. Ueland PM. Homocysteine Species as Components of Plasma Redox Thiol Status. Clin Chem 1995:41:340-342
- 2. Perry IJ. Refsum H. Morris RW. et al. Prospective Study of Serum Total Homocysteine Concentration and Risk of Stroke in Middle-aged British Men. The Lancet 1995;346:1395-1398
- 3. Refsum H. Smith AD. Ueland PM. et al. Facts and recommendations about total homocysteine determinations:an expert opinion. Clin Chem 2004:50(1):3-32
- 4. Nehler MR, Taylor LM Jr, Porter JM. Homocysteinemia as a Risk Factor for Atherosclerosis: A Review. Cardiovascular Pathol 1997;6:1-9
- 5. Mudd SH, Levy HL, Skovby F. Disorders of Transsulfuration. In: Scriver CR, Beaudet AL, Sly WS, et al., eds The Metabolic and Molecular Basis of Inherited Disease, New York: McGraw-Hill, 1995;1279-1327
- 6. Ueland PM, Refsum H, Stabler SP, et al. Total Homocysteine in Plasma or Serum: Methods and Clinical Applications. Clin Chem 1993;39:1764-1779