

ENZYMES

BY

SORACHIM

Peroxidase from Horseradish

PEO-302

SPECIFICATIONS

Product name	Donor:hydrogen-peroxidase oxidoreductase
EC	1.11.1.7
Appearance	Reddish brown amorphous powder lyophilized
Activity	Grade III, 180 Purpurogallin U/mg-solid or more (RZ \geq 2.0, salt free)
Contaminants	Phosphatase: < 0.001 %
Stability	Stable at -20 °C for at least 12 months
Molecular weight	approx. 40,000
Structure	Glycoprotein with one mole of protohaemin IX
Inhibitors	Cyanide, sulfide, fluoride, azide
Optimum pH	6.0 – 7.0
Optimum temperature	45 °C
pH Stability	pH 5.0 – 10.0 (25 °C, 20hr)
Thermal stability	below 50 °C (pH 6.0, 10min)

ENZYMES

BY

SORACHIM

Peroxidase from Horseradish

PEO-302

SPECIFICATIONS




Fig.1. Stability (Powder form)
(kept under dry conditions)




Fig.3. pH-Activity
(20°C, 20sec-reaction in 0.1M buffer solution: pH4.0-6.0, acetate; pH6.0-8.0, phosphate)




Fig.5. pH-Stability
(25°C, 20hr-treatment with 50mM buffer solution: pH3.5-6.0, acetate; pH6.0-8.0, phosphate; pH9.0-11.0 borate)




Fig.2. Stability (Powder form)
(kept under dry conditions)




Fig.4. Temperature activity
(20sec-reaction in 0.1M phosphate buffer, pH6.0)




Fig.6. Thermal stability
(10min-treatment with 50mM phosphate buffer, pH6.0)

E
N
Z
Y
M
E
S